FitSpine FX-2™ Inversion Table

with Comfort Cushion

For gravity-assisted stretching and decompression

FASTER, BETTER
RESULTS WITH
TEETER SUPPORT
& COACHING



ASSEMBLE FASTER

Free BILT app puts easy 3-D instructions at your fingertips. Search 'Teeter FX-2' in BILT to get started.



Find BILT in App Store



EASY ACCESS TO VIDEO

Watch the Getting Started DVD or access videos online for easy assembly guidance, how-to's and classes.



Scan or

Go to teeter.com/videos



FREE EXERCISE GUIDE

Register your warranty and gain access to our 30-Day Back Pain Relief Guide.



teeter.com/warranty



REACH YOUR GOALS

Use the free Teeter app to access guided sessions, product support, and more!





WELCOME TO THE TEETER FAMILY

Congratulations on your purchase of a Teeter FitSpine FX-2[™] Inversion Table! For the best experience, it is critical that you follow the assembly instructions, and read and fully understand the Owner's Manual attached to the equipment before inverting. Teeter Decompression Devices are multiple user, reusable medical devices for home use, intended to provide traction to the spine while stretching the para-spinal muscle and soft tissues. The devices provide non-powered traction and are meant for use by adults.

General Sector of 1 created Teeter so people could live healthier and more active lives. 99 -Roger Teeter

Teeter inversion tables are FDA-Registered as 510(k) medical devices. Indicated for:



- Back pain
- Muscle tension
- Degenerative disc disease
- Spinal degenerative joint disease
 Muscle spasm
- Spinal stenosis

- Herniated disc
- · Spinal curvature due to tight muscles
- Sciatica
- Facet syndrome

We've Got Your Back!



If you have any questions concerning assembly or if any parts are missing, DO NOT RETURN THE ITEM TO THE STORE OR CONTACT THE RETAILER. Our dedicated product service experts can help! Contact Teeter Customer Service at 800-847-0143, or via online forms or Live Chat at teeter.com.



Trust Teeter for unmatched quality and performance, backed by our industry-leading warranty coverage. To register your product warranty, go to **teeter.com/warranty**

TABLE OF CONTENTS

Assembly Instructions

Important Safety Instructions
Items for Assembly
A Guide to Your Inversion Table
Safety Warning Labels & Product Specifications 4
Assembly Steps5 – 13
Before Inverting
Warranty Terms & Registration

Owner's Manual (attached to the equipment)

Important Safety Instructions	
User Settings	2
Prepare to Invert	- 4
Inverting 4-	_ [
Storage & Maintenance	_
Get the Most Out of Your Teeter.	6

BEFORE YOU BEGIN: Review all steps before beginning assembly and read all precautions before using the inversion table. Carefully adhere to the Assembly Instructions and Owner's Manual to help ensure safety and product integrity.

IMPORTANT SAFETY INSTRUCTIONS READ ALL INSTRUCTIONS BEFORE USING THE INVERSION TABLE

AWARNING

FAILURE TO FOLLOW INSTRUCTIONS AND WARNINGS COULD RESULT IN SERIOUS INJURY OR DEATH.

To reduce the risk of injury:

- Read and understand all the instructions, review all other accompanying documents, and inspect the equipment before using the inversion table. It is your responsibility to familiarize yourself with the proper use of this equipment and the inherent risks of inversion if these instructions are not followed, such as falling on your head or neck, pinching, entrapment, equipment failure, or aggravating a pre-existing medical condition. It is the responsibility of the owner to ensure that all users of the product are fully informed about the proper use of the equipment and all safety precautions.
- **DO NOT** use until approved by a licensed physician. Inversion is contraindicated in any medical or health condition that may be made more severe by an elevation of blood pressure, intracranial pressure or mechanical stress of the inverted position, or that may impact your ability to operate the equipment. This may include injury or illness, but also the side effects of any drug or supplement (prescribed or over-the-counter). Specific conditions may include, but not be limited to:
 - Any condition, neurological or otherwise, which results in unexplained tingling, weakness or neuropathy, seizure, sleep disorder, lightheadedness, dizziness, disorientation, or fatigue, or impacts strength, mobility, alertness, or cognitive ability;
 - · Any brain condition, such as trauma, history of intracranial bleed, history or risk of TIA or stroke, or severe headaches;
 - Any condition of the heart or circulatory system, such as high blood pressure, hypertension, increased risk of stroke, or use of anticoagulants (including high doses of aspirin);
 - Any bone, skeletal or spinal cord condition or injury, such as significant spinal curvature, acutely swollen joints, osteoporosis, fractures, dislocations, medullary pins or surgically implanted orthopedic supports;
 - Any eye, ear, nasal or balance condition, such as trauma, history of retinal detachment, glaucoma, optic hypertension, chronic sinusitis, middle or inner ear disease, motion sickness, or vertigo;
 - · Any digestive or internal condition, such as severe acid reflux, hiatal or other hernia, gallbladder or kidney disease;
 - Any condition for which exercise is specifically directed, limited or prohibited by a physician, such as pregnancy, obesity, or recent surgery.
- **ALWAYS** be certain the Ankle Lock System is properly adjusted and fully engaged, and that your ankles are secure before using the equipment. HEAR, FEEL, SEE and TEST that the Ankle Lock System is snug, close-fitting and secure EVERY TIME you use the equipment
- · ALWAYS wear securely tied lace-up shoes with a flat sole, such as a normal tennis-style shoe.
- **DO NOT** wear any footwear that could interfere with securing the Ankle Lock System, such as shoes with thick soles, boots, high-tops or any shoe that extends above the anklebone.
- **DO NOT** use the inversion table until it is adjusted properly for your height and body weight. Improper settings can cause rapid inversion or make returning upright difficult. New users, and users who are physically or mentally compromised, will require the assistance of a spotter. Make sure the equipment is set to your unique user settings prior to each use.
- **DO NOT** sit up or raise head to return upright. Instead, bend knees and slide your body to the foot-end of the inversion table to change weight distribution. If locked out in full inversion, follow the instructions for releasing from the locked position before returning upright.
- **DO NOT** continue using the equipment if you feel pain or become light-headed or dizzy while inverting. Immediately return to the upright position for recovery and eventual dismount.
- **DO NOT** use if you are over 198 cm / 6 ft 6 in or over 136 kg (300 lb). Structural failure could occur or head / neck may impact the floor during inversion.
- **DO NOT** allow children to use this machine. Keep children, bystanders, and pets away from machine while in use. The inversion table is not intended for use by persons with reduced physical, sensory or mental capabilities, unless they are given supervision and instruction concerning use of the machine by a person responsible for their safety.
- DO NOT store the inversion table upright if children are present. Fold and lay the table on the floor. DO NOT store outdoors.
- **DO NOT** use aggressive movements, or use weights, elastic bands, any other exercise or stretching device or non-Teeter® attachments while on the inversion table. Use the inversion table only for its intended use as described in this manual.
- · DO NOT drop or insert any object into any opening. Keep body parts, hair, loose clothing and jewelry clear of all moving parts.
- · DO NOT use in any commercial, rental or institutional setting. This product is intended for indoor, home-use only.
- · DO NOT operate equipment while under the influence of drugs, alcohol, or medication that may cause drowsiness or disorientation.
- · ALWAYS inspect the equipment prior to use. Make sure all fasteners are secure.
- · ALWAYS replace defective components immediately and / or keep the equipment out of use until repair.
- · ALWAYS position equipment on a level surface and away from water or ledges that could lead to accidental immersion or falls.
- Refer to additional warning notices posted on the equipment. If a product label or Owner's Manual should become lost, damaged or illegible, contact Customer Service at 800-847-0143 for replacement.

ITEMS FOR ASSEMBLY

Parts are not shown to scale.

ITEM NO. ITEM NAME

Stretch-and-Grip™ A-Frame Base Assembly

E61100 A-Frame

F51007 Angle Tether pre-assembled to A-Frame

Handle Assembly

E61500 Stretch Assist[™] Handles (2)

Handle Assembly Hardware Kit HK1008

F50071 Allen Head Bolt (6)

Roller Hinge Assembly

F51064B Roller Hinges (2)

FitSpine™ Table Bed Assembly

EX1300 FitSpine™ Table Bed

FitSpine™ Table Bed Hardware Kit HK1010

EP1127 Bolt (2), H13009 Nut (2), EP1128 Washer (2)

ITEM NO. ITEM NAME

Main Shaft Assembly

EX1630 with EZ-Reach™ Ankle Lock System

Tools Provided for Assembly

IA1149 5mm Allen Wrench (1) EP1128A 6mm Allen Wrench (1) F51088 Open-Ended Wrench (1)

Product Support


112100 Owner's Manual pre-assembled to A-Frame

Included Accessories (Optional)

EX1350 Lumbar Bridge EX1450 Comfort Cushion EX1397 Attachment Nodes


Stretch-and-Grip™ A-Frame Base Assembly

with pre-assembled Angle Tether



FitSpine™ Table Bed Assembly

Use with Table Bed Hardware Kit (HK1010)







H13009

EP1128

Product Support



Owner's Manual

pre-assembled to A-Frame

1.12100

Handle Assembly

Use with Handle Assembly Hardware Kit (HK1008)



E61500







F50071

Main Shaft Assembly

EZ-Reach™ Ankle Lock System



Roller Hinge Assembly



F51064B

Tools Provided for Assembly



5mm Allen Wrench ΙΔ1149



6mm Allen Wrench EP1128A



Open-Ended Wrench

F51088

Included Accessories



Lumbar Bridge

EX1350



Comfort Cushion FX1450




Attachment Nodes

FX1397

A GUIDE TO YOUR INVERSION TABLE

Before reading further, study the drawing below to familiarize yourself with the important components of your Teeter Inversion Table.




SAFETY WARNING LABELS & PRODUCT SPECIFICATIONS

Important: Please review all labels and supporting materials before using your inversion table.

This drawing indicates the locations of the warning labels found on your product. If a label is missing, illegible or is removed, contact Teeter Customer Service to request a complimentary replacement label.

Note: Image and labels below not shown at actual size.



Assembled Non-Use Dimensions: $147.3 \times 73.1 \times 153.6 \text{ cm} / 58 \times 28.8 \times 60.5 \text{ in}$ Maximum In-Use Dimensions: $205.7 \times 73.1 \times 219.7 \text{ cm} / 81 \times 28.8 \times 86.5 \text{ in}$ Storage Dimensions: $50.8 \times 73.1 \times 167.6 \text{ cm} / 20 \times 28.8 \times 66 \text{ in}$

Weight (approx.): 28.4 kg / 62.6 lb

LET'S GET STARTED

Unpack and Prepare Your Workspace

- ·If possible, assemble the equipment at or near the space in which you intend to use it to avoid moving it later.
- · Unpack all parts and support materials. Set aside packing materials and clear your work area.
- ·Locate the Hardware Kits packaged with the manuals. They are labeled to correspond with the assembly process.


STEP 1

ASSEMBLE A-FRAME & HANDLES

Figure 1: On a level surface, position the A-Frame so that it is standing upright and the Stability Feet are on the ground.

Gently push down on the Spreader Arms to ensure they are fully open and in the "locked" position.

Look for temporary circular assembly assistance labels on the A-Frame. **RIGHT, LEFT, FRONT,** and **REAR** indicate your position while <u>using</u> the equipment, not facing it. These labels can be removed easily upon completion of assembly.




A - Spreader Arms B - Crossbar

Figure 1

Locate the Handle Assembly Hardware Kit (HK1008).

Figure 2: Determine the left or right handles, marked with an embossed L/R on the inside of the black plastic part of each handle.

Align the black plastic part of the corresponding handle (left / right) over the outside edge of the Hinge Plate on the A-Frame.



Repeat with other handle. Tighten all fasteners with the 5mm Allen Wrench provided, being careful not to over-tighten.




Figure 2




Figure 3

ASSEMBLE TABLE BED




Figure 4

Locate the Table Bed Assembly Hardware Kit (HK1010).

Figure 4: Place the Table Bed face down on the floor and push down on the Support Beam, so the two holes align evenly with the holes at the base of the Upper Support Arms. You may have to exert extra pressure to ensure that the Support Beam slides over the rubber spacers.




Figure 5

Figure 5: Insert the two Bolts into the open holes.

Thread a Washer and Nut onto each Bolt and hand-tighten.

Using the 6mm Allen Wrench to steady the bolts, tighten the Nuts onto the Bolts with the 10/13mm Open-Ended Wrench.

ASSEMBLE ROLLER HINGES TO TABLE BED

Figure 6: Familiarize yourself with the 3-Hole Roller Hinge and Cam Lock terms.

AWARNING

NEVER disassemble the Roller Hinge Pivot Pin.
ALWAYS insert the 3-Hole Roller Hinge (with the Pivot Pin on top and facing out) in the same direction as the arrow label located inside of the Cam Lock for proper assembly.




Figure 7a: On one side of the Table Bed, lift and hold the Cam Lock up all the way to unlock.

In your other hand, hold one Roller Hinge near the Pivot Pin. With the Pivot Pin facing out (away from the Table Bed), slide the bottom of the Roller Hinge between the Cam Lock and the Bracket in the same direction as the arrow label located inside of the Cam Lock.


TIP: Make sure that the Cam Lock is completely open when inserting the Roller Hinge, otherwise assembly will be more difficult.



Repeat on other side. Make sure the Roller Hinges are in the same hole setting on both sides.



NOTE: Refer to the Owner's Manual for an explanation of the hole settings. If you are unsure, use Setting C to start.






Figure 7a




Figure 8




Figure 9

ASSEMBLE TABLE BED TO A-FRAME




Figure 10

Figure 10: Face the front of the A-Frame where the Crossbar is located. Grasp both Roller Hinges, right above the Cam Lock, and lift the Table Bed. Allow the top of the Table Bed to rotate toward the floor, so that the back of the Table Bed is now facing you and the top of the Table Bed is in front of the Crossbar.




Figure 11




Figure 11a




Figure 11b

Figure 11: Lower each Roller Hinge Pivot Pin into the A-Frame hinge plates, one side at a time. The Self-Locking Hooks will open to allow the Pivot Pin into the Hinge Plate slot, then automatically snap closed over the Pivot Pin.

TIP: You may need to push outward on the Hinge Plate in order for the second Pivot Pin to lock in place.

Figure 11a & 11b: Make sure that each Pivot Pin is seated at the base of the slot in the Hinge Plates, and that the Self-Locking Hooks have closed over both Pivot Pins.

AWARNING

Failure of the Self-Locking Hooks to close over both Roller Hinge Pivot Pins is an indication of improper assembly and if not corrected could result in serious injury or death!

Figure 12: Rotate the Table Bed so that it is facing up. Ensure that it rotates smoothly. See also Image A on Page

10 to ensure correct assembly.




Figure 12

ASSEMBLE MAIN SHAFT TO TABLE BED

Figure 13: Stand on the LEFT side of the A-Frame, holding the Main Shaft with the height markings facing up. Begin to slide the end of the Main Shaft into the Main Shaft Housing at the base of the Table Bed.

Figure 13a: With your right hand, pull out the Height-Selector Locking Pin and slide the Main Shaft in further. Release in the desired height setting. Refer to the Owner's Manual for more information on selecting your height setting.




Figure 13

Figure 14: The Main Shaft MUST REST against the Crossbar of the A-Frame.

IMPORTANT: The Crossbar prevents the Table Bed from rotating forward when the user steps on the Ankle Comfort Dial. If the Main Shaft does not rest on the Crossbar as shown here, then the Table Bed has been assembled backwards onto the A-Frame. This must be corrected before use.





Figure 14

Figure 15: Test the inversion table by hand for smooth and steady rotation and ensure that all fasteners are secure.



MISASSEMBLY CHECK

AWARNING

If your Teeter Inversion Table looks like either of these images, your inversion table has been misassembled and is unfit for use. Improper assembly could result in serious injury or death!

Image A

Go back to Step 3 for instruction.

Demonstrates that the Roller Hinges have been assembled upside down into the Table Bed and must be corrected.





Image B

Go back to Step 4 for instruction.

Demonstrates that the Table Bed has been assembled into the A-Frame backwards so the Main Shaft is not resting on the Crossbar and must be corrected.



STEP 6 ATTACH ANGLE TETHER

Attach Angle Tether

The tether will come pre-assembled to the A-Frame.

Figure 16: Unfold the adjustable tether and clip it to the U-Bar on the underside of the Table Bed.

Slide the buckle to lengthen or shorten the strap depending on your desired maximum angle of inversion.




Figure 16

ATTACH COMFORT CUSHION (OPTIONAL)

The Comfort Cushion provides a plush, supportive surface for a premium inversion experience without inhibiting body slide needed for decompression. The included fastening nodes secure the Comfort Cushion to the Table Bed while still allowing it to move within the tracks of the FitSpine Bed.



Figure 17

AWARNING

CHOKING HAZARD: Small parts - keep out of reach of children.

Figure 17: Locate the fastening nodes included with the Comfort Cushion and unscrew the backing from each.

Place the cushion onto the Table Bed aligned with the designated holes.



Figure 18



Figure 18a



Figure 18b

Figure 18: Insert a node backing from the back of the Table Bed through the Table Bed slot and hole in the cushion. Hold the backing perpendicular to the slot.

Figure 18a & 18b: With the other hand, fasten the node top onto the backing threads visible from the front of the Table Bed. Twist the node clockwise to secure it into place in one of two ways:

- 1. Static: If you do not want the Comfort Cushion to move with you as you invert, tighten the node until it is snug. Do not over-tighten.
- 2. Responsive: If you want the Comfort Cushion to slide with you as you invert to encourage body slide at shallower angles, tighten the node just tight enough to keep it in place, but loose enough to move within the slot.

Repeat the process above with the three (3) remaining fastening nodes to finish securing the Comfort Cushion to the Table Bed.

Be sure to replace the backings onto the nodes when the Comfort Cushion is removed.

Cleaning & Maintenance

Spot clean using a soft, slightly damp sponge or cloth. Do not use brushes or abrasive cleaners.

ATTACH LUMBAR BRIDGE (OPTIONAL)

The Lumbar Bridge provides even deeper decompression and improved alignment benefits. Personalize the intensity and target zone of the Lumbar Bridge by adjusting the height and position of the arch within the slots on the Table Bed so it fits comfortably at the small of your back. The Lumbar Bridge can be used with or without the Comfort Cushion.

Stabilize the inversion table to prevent rotation during assembly.

Position the bridge with the Teeter logo at the base and facing towards you.

Figure 19: First, insert the bottom two notches of the bridge into the desired horizontal slots on the lower portion of the Table Bed.



Figure 19

Figure 20: Bend the bridge to the arch height you desire and insert the one top notch into the corresponding slot on the upper portion of the Table Bed.

Modify the positioning and arch height as needed.

NOTE: Long-term storage of the Lumbar Bridge in a high arch setting may result in distortion of the shape and ability of the bridge to move to lower arch settings. Store in its flat position when not in use.



Figure 20

Use Instructions

Use of the Lumbar Bridge may impact the rotation of the table. Set the Angle Tether to a moderate inversion angle and use a spotter until you are able to find the correct balance settings and are comfortable with the operation of the table.

Start in a lower level arch setting and work your way up. If you feel any discomfort, lower the setting until you are comfortable or discontinue use.

Figure 21: The Lumbar Bridge works between intermediate to moderate inversion. Full inversion moves the body away from the table bed and will reduce the effectiveness of the bridge.



Figure 21

BEFORE INVERTING




Figure 22




Figure 22a

Ensure Owner's Manual is Attached

The Owner's Manual contains important information on how to use your Teeter Inversion Table, including how to personalize the user settings, properly secure and release the Ankle Lock System, and test and adjust the rotation control.

Figure 22: If the Owner's Manual is not already attached, secure the chain to the A-Frame through the designated hole in the Hinge Plate. Allow the Owner's Manual to hang freely on the outside of the A-Frame Spreader Arms so it doesn't interfere with the rotation of the Table Bed.

IMPORTANT: Once attached to the A-Frame, DO NOT remove the Owner's Manual. It should remain permanently attached to your inversion table to serve as a reference for all users in regard to proper adjustment and use of the equipment.

AWARNING

Read the Owner's Manual thoroughly before using your Teeter Inversion Table. Improper settings could result in serious injury or death!

Watch the Getting Started Videos

The Getting Started Video Portal is a helpful supplement to the Owner's Manual, with easy-to-follow instructions on user settings, how to invert, storage and maintenance, and even stretching and exercises you can do with your Teeter. Access via teeter.com/videos and search Model 'FX-2' now!

The Teeter warranty set forth below and on Teeter's website applies to US and Canadian customers only. For international customers, please consult your local distributor for warranty information which will vary depending on country.



principal place of business.

During the period starting with the day of retail purchase and continuing for five (5) years, Teeter extends to the owner a repair and replacement warranty against manufacturing defects in materials, workmanship, fabrics and padding. Teeter will repair or replace any such defect and will pay the costs of all parts, labor and transportation. If a repair or replacement is not commercially practical or cannot timely be made, then Teeter will, at the original Purchaser's option, replace with a comparable product or refund the purchase price.

Handling and transportation costs related to product warranty service only are covered by this warranty. This warranty does not cover damage resulting from improper handling, assembly, or installation, repairs made by others, accident, misuse, or abuse. Under no circumstances shall Teeter, or any other party involved in the sale of this product, have any liability for incidental or consequential damage arising from breach of an express or implied warranty on any Teeter product.

EXCEPT AS SET FORTH ABOVE, NO WARRANTY IS GIVEN WITH RESPECT TO ANY TEETER PRODUCT, AND ALL EXPRESS WARRANTIES ARE DISCLAIMED. This warranty shall be governed by the laws of the State of Washington, USA. To the extent this warranty is found not to be enforceable, it shall be deemed revised to the extent necessary to make it enforceable. This warranty and any controversy or claim arising out of this warranty or its interpretation shall be governed by the laws of the State of Washington, USA. Any controversy or claim arising out of or relating to this warranty, its interpretation, or any alleged breach thereof, which cannot be amicably settled between Teeter and the owner within sixty (60) days of written notice by the aggrieved party to the other, shall be finally settled by arbitration submitted to three (3) arbitrators selected from the panels of the arbitrators of the American Arbitration Association located closest to Teeter's

Some states do not allow the exclusion of incidental or consequential damage from a warranty, so the above limitation or exclusion may not apply to you. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

AUSTRALIA: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

HOW TO SUBMIT YOUR REGISTRATION:

Step 1

Fill out this information for your own records.

Step 2

Go online to teeter.com/warranty to register your warranty.

Date of Purchase

Product & Model

Dealer Name

Serial No. (located on the back of the Table Bed)



FREE DIGITAL GUIDE AVAILABLE!

Maximize your results with your Teeter Inversion Table

Go to teeter.com and register your warranty to gain **FREE** access to this downloadable exercise guide.

If you are unable to go online, you can request a warranty card to be mailed to you by calling Customer Service at 800-847-0143.

Please DO NOT mail this to Teeter.







Distributed by:
Costco Wholesale Corporation
P.O. Box 34535
Seattle, WA 98124-1535
USA
1-800-774-2678
www.costco.com

Costco Wholesale Canada Ltd.*
415 W. Hunt Club Road
Ottawa, Ontario
K2E 1C5, Canada
1-800-463-3783
www.costco.ca

* faisant affaire au Québec sous le nom les Entrepôts Costco

Importado por:
Importadora Primex S.A. de C.V.
Blvd. Magnocentro No. 4
San Fernando La Herradura
Huixquilucan, Estado de México
C.P. 52765
RFC: IPR-930907-S70
(55)-5246-5500
www.costco.com.mx

Costco Wholesale Iceland ehf.

Kauptún 3-7, 210 Gardabaer
Iceland
Www.costco is

Costco France
1 avenue de Bréhat
91140 Villebon-sur-Yvette
France

Costco Wholesale UK Ltd / Costco Online UK Ltd Hartspring Lane Watford, Herts WD25 8JS United Kingdom 01923 213113 www.costco.co.uk

Costco Wholesale Japan Ltd. 3-1-4 Ikegami-Shincho Kawasaki-ku, Kawasaki-shi, Kanagawa 210-0832 Japan 0570-032600

www.costco.co.jp

Costco Wholesale Spain S.L.U.
Polígono Empresarial Los Gavilanes
C/ Agustín de Betancourt, 17
28906 Getafe (Madrid) España
NIF: B86509460
900 111 155
www.costco.es

Costco Wholesale Korea, Ltd. 40, Iljik-ro Gwangmyeong-si Gyeonggi-do, 14347, Korea 1899-9900 www.costco.co.kr Costco President Taiwan, Inc. No. 656 Chung-Hwa 5th Road Kaohsiung, Taiwan Company Tax ID: 96972798 449-9909 or 02-449-9909 (if cellphone) www.costco.com.tw Shanghai Minhang Costco Trading Co., Ltd No. 235, ZhuJian Road Minhang District, Shanghai China 201106 +86-21-6257-7065

01 80 45 01 10

www.costco.fr

Costco Wholesale Australia Pty Ltd

17-21 Parramatta Road

Lidcombe NSW 2141

Australia

www.costco.com.au

ATEETER

If you have any trouble assembling the equipment, or questions about its use, please contact customer service.

USA: 800-847-0143 or info@teeter.com (English language) International: info@teeterintl.com

Hours of Operation (Monday - Friday 8am - 4pm PST)

USA: Teeter 9713 233rd Avenue East Ste A Bonney Lake, WA 98391 Toll Free: 800-847-0143 Fax: 800-847-0188 teeter.com | info@teeter.com



Any modification to this device will void the UL Listing



Medical Device Safety Service GmbH Schiffgraben 41 30175 Hannover Germany Tel. +49 511 62628630



Teeter Decompression Devices are multiple user, reusable devices for home use, intended to provide traction to the spine while stretching the para-spinal muscles and soft tissues.

The devices provide non-powered traction and are meant for use by adults.

Use of the Teeter Decompression Devices is indicated for the following conditions: back pain, muscle tension, degenerative disc disease, spinal degenerative joint disease, spinal stenosis, herniated disc, spinal curvature due to tight muscles, sciatica, muscle spasm, and facet syndrome.